

# HORIZONTAL AUTOCLAVE AUT22-0146SDA



# HORIZONTAL AUTOCLAVE AUT22-0146SDA

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## AUT22-0146SDA HORIZONTAL AUTOCLAVE



Product Image Coming Soon

### SPECIFICATIONS

|                         |                   |
|-------------------------|-------------------|
| Model                   | AUT22-0146SDA     |
| Capacity                | 146 L             |
| Chamber Size (Wx H x L) | 420x420x830 mm    |
| External Size (WxHxL)   | 1290x1680x1100 mm |
| Net Weight              | 610 Kg            |
| Power Consumption       | 18 KW             |
| Power Supply            | 380 V 50 Hz       |



**Labtare USA**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)