

HORIZONTAL AUTOCLAVE AUT22-0665SDB




HORIZONTAL AUTOCLAVE AUT22-0665SDB

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

AUT22-0665SDB HORIZONTAL AUTOCLAVE



Product Image Coming Soon

SPECIFICATIONS

Model	AUT22-0665SDB
Capacity	655 L
Chamber Size (Wx H x L)	610x910x1180 mm
External Size (WxHxL)	1310x2070x1420 mm
Net Weight	1370 Kg
Power Consumption	46 KW
Power Supply	380 V 50 Hz



Labtare USA

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@labtare.com | Website: labtare.com