

MICROPLATE CENTRIFUGE CEN31-0250MP



MICROPLATE CENTRIFUGE CEN31-0250MP

These laboratory multi purpose Medical Centrifuges are ideal solution for high-end research experiment such as PRP/PPP, cell washing, HCT, blood type test etc. in the biology, radio-immunity, medical and other research fields. These versatile centrifuges have precise temperature control with highly efficient double loop control design which allows fast heating and cooling.

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover. Also known as Benchtop Centrifuge, Laboratory Centrifuge.

CEN31-0250MP MICROPLATE CENTRIFUGE

Applicable to all capacity PCR plates including 96-hole/384-hole microplate



SPECIFICATIONS

Model	CEN31-0250MP
Maximum Capacity (No of tubes x Vol.)	384-well PCR microliter plate
Maximum Speed RPM	2500
Maximum RCF	500xg
Ambient Temperature	5°C~40°C
Relative Humidity	≤85%
Atmospheric Pressure	860hPa-1060hPa
Overall Dimension	222x220x190 mm
Noise Level	≤ 60 dB
Weight (Net/Gross)	2.3 kg
Power Supply	110V 50/60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Capacity
LS55168	Microplate Rotor	384-well / 96-well / Small capacity microliter plates / PCR microliter plates with or without skirt	384-well Microplate Rotor



Labtare USA

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@labtare.com | Website: labtare.com