

HEMATOLOGY ANALYZER LTANA10-2



HEMATOLOGY ANALYZER LTANA10-2

Also known as Blood Analyzer, Laboratory Hematology Analyzer.

LTANA10-2 HEMATOLOGY ANALYZER

High Efficiency and convenience: (1) 60 samples per hour, minimize your workload. (2) 1 click for analysis with 20 parameters+3 histograms.

Complete diagnostic solution: (1) Two reagents only(500mL LYSE, 20 L DILUENT), cost-effective. (2) International brand control and calibrator provided, high quality.

Powerful data management: (1) Store 100,000 results, large memory. (2) Easy to check history results, humanized design. (3) Support LIS, easy data transmission.

Intelligent design: (1) Intelligent software menu for hardware diagnosis. (2) Advanced high voltage and flush to remove blockage automatically.



SPECIFICATIONS

Model	LTANA10-2
Interface	USB port, RS232
Test Principle	Impedance for cell counting and cyanide-free method for HGB
Clog Clear	High voltage and high pressure flush
Parameter	WBC, Neu#, Lym#, Mid#, Neu%,Lym%, Mid%, RBC, HGB, HCT, MCV, MCH, MCHC,RDW-SD, RDW-CV, PLT, MPV, PDW, PCT, P-LCR
Histograms	WBC, RBC, PLT
Sample mode	Venous mode: 10pL, venous blood Capillary mode: 10pL, capillary blood Prediluted mode: 20pL capillary blood
Language	English, Spanish, Portuguese, Italian, etc.
QC Mode	3 level QC, L-J graph
Calibration Mode	Auto mode and Manual mode
Reagent	Lyse(500ml) and Diluent(20L) Reagent required for single test: Diluent: 26mL, Lyse: 0.35mL
Memory	100,000 sample results
Display	Touch screen LCD
Printer	Thermal printer, supports external printer
Dimension	514x594x645 mm
Weight	24.7 kg
Power Supply	100~240 AC, 50/60±1 Hz

ACCESSORIES

Accessory Code	Name	Unit	1
1801008020	Power Cable (1.8 m)	pc	1
1801008021	Ground Cable	pc	1
1801008022	Serial Port Cable	pc	1
1801008023	Mouse	pc	1
1801008024	Keyboard	pc	1

1801008025	Touch Pen	pc	1
1801008026	Side Door Key	pc	1
1801008027	Hex Wrench (Φ 2.5)	set	1
1801008028	Service Pack	set	1
1801008029	Thermal Printer Paper (Width: 80 mm)	roll	1
1801008030	Diluent (20 L)	bottle	1
1801008031	Lyse Solution (500 mL)	bottle	1
1801008032	PB Cleanser (50 mL)	bottle	1
1801008033	E-Z Cleanser (50 mL)	bottle	1
1801008034	Diluent Cap Component	set	1
1801008035	Lyse Solution Cap Component	set	1
1801008036	Waste Bottle Cap Component	set	1
1801008037	Waste Bottle	pc	1

OPTIONAL ACCESSORIES

Accessory Code	Name	Unit
1801008008	Thermal Printer Paper (Width: 80 mm)	roll
1801008009	PB Cleanser (50 mL)	bottle
1801008010	E-Z Cleanser (50 mL)	bottle
1801008011	Service Manual	pc
1801008012	Diluent (20 L)	bottle
1801008013	Lyse solution (500 mL)	bottle

Accessory Code	Name	Unit
1801008014	CBC-3D Blood QC High Level (3.0 mL) (R&D)	tube
1801008015	CBC-3D Blood QC Normal Level (3.0 mL) (R&D)	tube
1801008016	CBC-3D Blood QC Low Level (3.0 mL) (R&D)	tube
1801008017	Lyse solution (1 L)	bottle
1801008018	CBC-CAL PLUS Calibrator (3.0 mL) (R&D)	tube
1801008019	Barcode Scanner Assembly	pc



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com