

# HALOGEN MOISTURE ANALYZER LTANA12-5



# HALOGEN MOISTURE ANALYZER LTANA12-5

Also known as Laboratory Moisture Analyzer, Moisture Balance, Moisture Meter, Halogen Moisture analyzer, Laboratory moisture Meter, Laboratory Halogen Moisture analyzer.

## LTANA12-5 HALOGEN MOISTURE ANALYZER

Halogen Light Heating

Stainless Steel Chamber

Temperature and time can be set

The percentage of moisture content

The percentage of dry residual

The super high resolution back window matrix liquid crystal display (LCD), making it easier to operate in the dark place, and the user will have a more comfortable vision.

Stores historical setting



## SPECIFICATIONS

Model	LTANA12-5
Capacity	110g/5mg
Temperature	40°C-199°C
Operating Temperature Range	5°C-35°C
Temperature Set	1°C
Moisture Range	0.00%-100.0%
Moisture Readability	0.001
Dry Range	100.00%-0.00%
Dry Readability	0.001
Pan Size	Φ90 mm
Readability	5 mg
Heat Source	HALOGEN LAMP-1
ATRO	100%....999%
ATRO 2	0%....999%
Sensor	LOADCELL
Temperature Sensor	PT-100
Packaging Size	490x360x350 mm
Number of Storage	15
Weight	7.5 kg



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)