

# AUTOSAMPLER FOR TOTAL ORGANIC CARBON TOC ANALYZER LTANA13-7



# AUTOSAMPLER FOR TOTAL ORGANIC CARBON TOC ANALYZER LTANA13-7

## LTANA13-7 AUTOSAMPLER FOR TOTAL ORGANIC CARBON TOC ANALYZER

Compact design, only occupy little space.

20 positions, meet most of experiment needs.

Easy to operate, no need special operation training.

Equipped with USB port for PC connection.

Equipped with multi-sensors for real-time monitoring which ensures precision positioning and operation safety.

No need special calibration or maintenance after long time use.



### SPECIFICATIONS

Model	LTANA13-7
Sample Bottle Volume	60ml
Maximum Number	19
Injection Arm Stroke	85 mm
Ambient Temperature	0-40 °C
Relative Humidity	≤85%
Rated Power	120W
Power Adapter	AC 100-240V,50/60HZ
Dimension	265x373x395 mm

# LABTARE

**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)