

MILK ANALYZER LTANA14-1



MILK ANALYZER LTANA14-1

Also known as Milk Testing Instrument, Ultrasonic Milk Analyzer..

LTANA14-1 MILK ANALYZER



SPECIFICATIONS

Model	LTANA14-1
Measuring Range	
Fat	0.01%-25%
SNF	3%-15%
Density	1015-1040 kg/m ³
Proteins	2%-7%
Lactose	0.01%-6%
Water Content	0%-70%
Temperature	1°C-40°C
Freezing Point	-0.4°C- -0.7°C
Salt	0.4%- 1.5%
Accuracy	
Fat	± 0.10%
SNF	± 0.15%
Density	± 0.3 kg/m ³
Proteins	± 0.15%
Lactose	± 0.20%
Water Content	± 3.0%
Temperature	±1°C
Freezing Point	±0.001°C
Salts	± 0.05%
Interface	RS232
Standard Accessory	Printer
Display	LCD Display – 4 lines x 16 characters
Power	100-240 V , 50/60 Hz
Gross Weight	6.5 kg
Package Size (WxDXH)	30x30x45 cm



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com