

FLOOR STANDING HORIZONTAL AUTOCLAVE LTAUT14-2



FLOOR STANDING HORIZONTAL AUTOCLAVE

LTAUT14-2

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

LTAUT14-2 FLOOR STANDING HORIZONTAL AUTOCLAVE

High Quality Stainless steel chamber with drying function

Safety features for over temperature and pressure

Automatic opening of safety valve if pressure exceeds 0.165 MPa

Door remains locked if pressure exceeds 0.027MPa

Over temperature, over pressure, auto protection device

Automatic power cut off and water cut off in case of water shortage



SPECIFICATIONS

Model	LTAUT14-2
Capacity	200 L
Adjustment of Temperature	105°C - 134°C
Sterilization Temperature	134°C
Sterilization Time	0-60 min
Heat Average	$\leq \pm 1^{\circ}\text{C}$
Working Pressure	0.22 MPa
Time of Drying	0-60 min
Chamber Dimension	515x1000 mm
Overall Dimension (LxWxH)	1400x670x1650 mm
Packing Size	1560 x 820 x 1850 mm
Weight	260/350 kg
Power	9000 W
Power Supply	380 V, 50Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com