

# FLOOR STANDING HORIZONTAL AUTOCLAVE LTAUT19-2



# FLOOR STANDING HORIZONTAL AUTOCLAVE

## LTAUT19-2

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## LTAUT19-2 FLOOR STANDING HORIZONTAL AUTOCLAVE

Sterilizing programs for fabric, instrument, rubber and liquid

Touch screen display of temperature, pressure, time, running status and failure warning

Post-drying function makes the residual humidity of sterilized articles less than 3%

The positive pressure air-exhaust system to exhaust the cold air from chamber

Equipped with printer to print date, time and process parameters

Build-in efficient steam generator, microcomputer temperature control

Perfect four-level authority management system to avoid misoperation

Double door optional



### SPECIFICATIONS

|                         |                   |
|-------------------------|-------------------|
| Model                   | LTAUT19-2         |
| Capacity                | 135 L             |
| Chamber Size            | Φ500x710 mm       |
| Chamber Material        | SUS304            |
| Max Working Pressure    | 0.23 Mpa          |
| Vacuum Limit            | -0.08 Mpa         |
| Working Temperature.    | 115-138°C         |
| Temp Accuracy           | 10.1°C            |
| Consumption             | 11.5 KW           |
| Power Supply            | 380 V,50 Hz       |
| External Size(Lx Wx H)  | 1110x750x1742 mm  |
| Package Size(L x W x H) | 1370x1000x2060 mm |
| Gross Weight            | 470 kg            |



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)