

FLOOR STANDING HORIZONTAL AUTOCLAVE LTAUT19-3



FLOOR STANDING HORIZONTAL AUTOCLAVE

LTAUT19-3

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

LTAUT19-3 FLOOR STANDING HORIZONTAL AUTOCLAVE

Sterilizing programs for fabric, instrument, rubber and liquid

Touch screen display of temperature, pressure, time, running status and failure warning

Post-drying function makes the residual humidity of sterilized articles less than 3%

The positive pressure air-exhaust system to exhaust the cold air from chamber

Equipped with printer to print date, time and process parameters

Build-in efficient steam generator, microcomputer temperature control

Perfect four-level authority management system to avoid misoperation

Double door optional



SPECIFICATIONS

| | |
|-------------------------|-------------------|
| Model | LTAUT19-3 |
| Capacity | 185 L |
| Chamber Size | Φ500x950 mm |
| Chamber Material | SUS304 |
| Max Working Pressure | 0.23 Mpa |
| Vacuum Limit | -0.08 Mpa |
| Working Temperature. | 115-138°C |
| Temp Accuracy | 10.1°C |
| Consumption | 11.5 KW |
| Power Supply | 380 V,50 Hz |
| External Size(Lx Wx H) | 1350x750x 1742 mm |
| Package Size(L x W x H) | 1510x890x2060 mm |
| Gross Weight | 540 kg |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com