

FLOOR STANDING HORIZONTAL AUTOCLAVE LTAUT19-39



FLOOR STANDING HORIZONTAL AUTOCLAVE

LTAUT19-39

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

LTAUT19-39 FLOOR STANDING HORIZONTAL AUTOCLAVE



SPECIFICATIONS

Model	LTAUT19-39
Capacity	650 L
Chamber Size (Wx H x L)	672x672x1500 mm
External Size (WxHxL)	1314x1260x1980 mm
Max Designed Pressure	Inner chamber - 0.1/0.3 Mpa, Jacket 0.3MPa
Working Pressure	0.23 Mpa
Built - in Steam Generator	0.8 Mpa
Designed Temperature	150 °C
Maximum Working Temperature	134 °C
Temperature Control Precision	± 0.2 °C
Consumption	56kVA
Standard Accessories	Trolley, Printer
Altitude	≤ 2000m; could be customised @ 2000-3000m
Power Supply	380 V 50 Hz
Package Size	Main Instrument: 1424x1460x2295; Trolley: 1860x640x1245
Gross Weight	Main Instrument: 1440; Trolley: 130



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com