

# FLOOR STANDING HORIZONTAL AUTOCLAVE LTAUT19-4



# FLOOR STANDING HORIZONTAL AUTOCLAVE

## LTAUT19-4

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## LTAUT19-4 FLOOR STANDING HORIZONTAL AUTOCLAVE

Sterilizing programs for fabric, instrument, rubber and liquid

Touch screen display of temperature, pressure, time, running status and failure warning

Post-drying function makes the residual humidity of sterilized articles less than 3%

The positive pressure air-exhaust system to exhaust the cold air from chamber

Equipped with printer to print date, time and process parameters

Build-in efficient steam generator, microcomputer temperature control

Perfect four-level authority management system to avoid misoperation

Double door optional



### SPECIFICATIONS

Model	LTAUT19-4
Capacity	300 L
Chamber Size	Φ632x1000 mm
Chamber Material	SUS304
Max Working Pressure	0.23 Mpa
Vacuum Limit	-0.08 Mpa
Working Temperature.	115-138°C
Temp Accuracy	10.1°C
Consumption	18 KW
Power Supply	380 V,50 Hz
External Size(Lx Wx H)	1400x890x1780mm
Package Size(L x W x H)	1700x1050x2090 mm
Gross Weight	830 kg



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)