

BENCHTOP AUTOCLAVE LTAUT22-1



BENCHTOP AUTOCLAVE LTAUT22-1

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Medical Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

LTAUT22-1 BENCHTOP AUTOCLAVE



Product Image Coming Soon

SPECIFICATIONS

Model	LTAUT22-1
Capacity	3 L
Sterilization Pressure	0.22 Mp
Sterilization Temperature	134 °C
Sterilization Time	4 min
Sterilization Cycle time	11 min
Required water	Distilled water or purified water
Product Size	300X190X350 mm
Outer Dimension	370x260x460 mm
Gross Weight	7.5 Kg
Power	0.5 KW
Voltage	AC220/110V+10%



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com