

BENCHTOP AUTOCLAVE LTAUT22-4



BENCHTOP AUTOCLAVE LTAUT22-4

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Medical Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

LTAUT22-4 BENCHTOP AUTOCLAVE



Product Image Coming Soon

SPECIFICATIONS

Model	LTAUT22-4
Capacity	23 L
Product Size	450x550x350 mm
Chamber Dimension	φ249x450 mm
Package Size	670x550x500 mm
Gross Weight	33 Kg
Net Weight	26 Kg
Power	800 W
Voltage	AC220/110V+10%, 50/60 Hz+ 1Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com