

BENCHTOP AUTOCLAVE LTAUT24-1



BENCHTOP AUTOCLAVE LTAUT24-1

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

LTAUT24-1 BENCHTOP AUTOCLAVE

Pre and Post Vacuum Sterilizer

Dual core control system

11 Preset Sterilization Cycles

Maximum negative pressure of -0.80 bar

Preheated chamber keeps the chamber warm and ready to use

LCD display

High density bacteria filter

User-defined cycles and cleaning programs

Drying time: 1-60 min

Printer: Thermal/Stylus

Interface: RS232

ISO & CE 1366 Certified

Vacuum Helix and Bowie & Dick Test passed

Powerful steam generator and vacuum pump

First class safety system with Electrical lock

Multi languages available: English, Italian, German, Chinese

Separate and set up and service menu

Password protection

Safety features: Over heat control, over pressure control, safety door locking system, low water level indicator



SPECIFICATIONS

Model	LTAUT24-1
Capacity	8 L
Sterilization Temperature	121°C, 134°C
Sterilization Time	4 min - 60 min
Working Pressure	0.21 - 0.23 MPa
Water Tank	2 x 4.5L Internal Stainless Steel water tanks
USB output	yes
Lock system of door	Mechanic + Electrical lockers
Chamber Dimension	170x320 mm
Overall Dimension (LxWxH)	535x365x320 mm
Packing Size	780 x 500 x 430 mm
Chamber Material	Stainless Steel SUS 304
Weight	59/70 kg
Power	2400 W
Power Supply	220 V, 50Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Unit
1900808008	Ring	1
1900808009	Inner temperature sensor	1
1900808010	pt1000	1
1900808011	pt1000 (packing with plastic)	1
1900808012	Heating circle	1
1900808013	Steam generator cover	1Set
1900808014	Steam generator stand	1Set
1900808015	Steam generator spring type gasket	1
1900808016	Heating Rod	3
1900808017	Water Pump	1
1900808018	Vacuum Pump	1
1900808019	Fuse Wire	2
1900808020	Normally closed solenoid valve	3

Accessory Code	Name	Unit
1900808021	Three-way solenoid valve	2
1900808022	Germ-tight filter	1
1900808023	Voltage transformer	1
1900808024	switch	1
1900808025	Bleeder valve	2
1900808026	Fan	1
1900808027	Printer	1
1900808028	Control panel	1set
1900808029	LED	1set
1900808030	Big Pipe	1Kg
1900808031	Small Pipe	1kg
1900808032	Liquid sensor	3
1900808033	Electromagnet	1



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com