

# FLOOR STANDING VERTICAL AUTOCLAVE LTAUT30-6



# FLOOR STANDING VERTICAL AUTOCLAVE

## LTAUT30-6

Used in Food industry, Chemical industry, Research, Laboratory, Pharmaceutical.

Also known as Vertical Pressure Sterilizer, Vertical Pressure Autoclave, Vertical Sterilizer, Laboratory Vertical Autoclave,.

### LTAUT30-6 FLOOR STANDING VERTICAL AUTOCLAVE

Microcomputer control system

Unique anti scald design chamber for corrosion resistance

Stainless steel chamber

Five sterilization modes

User friendly control panel

Multi layered pressure protection

Fully automatic internal discharge control embedded with a steam trapping bottle

Automatic checking of locking condition of cover

Automatic sterilization according to preset time

Safety device: Self-induction pressure interlocking device, over temperature protection system, dry scorch protection system, overpressure protection, safety valve, over current and short circuit protection, cover inspection system, leakage protection device, anti-scald safety protection and automatic troubleshooting system



Product Image Coming Soon

### SPECIFICATIONS

Model	LTAUT30-6
Capacity	54 L
Function	Fully automatic, No flow chart
Temperature Range	5°C-138°C
Sterilization Temperature	105-135°C
Warming Temperature	45-60°C
Warming Time	1-9999 min
Melting Temperature	60-100
Melting Time	1-300 min
Humidity Range	10%-85% RH
Working Pressure	0-0.6 MPa
Chamber Dimension	320x730 mm
Overall Dimension (LxWxH)	460x542x1070 mm
Weight	80/120 kg
Power	2900 W
Power Supply	220 V, 60Hz

### ACCESSORIES

Accessory Code	Name
1903008008	Stainless steel baskets
1903008009	Drainage hose

1903008010	Exhaust bottle
------------	----------------

## OPTIONAL ACCESSORIES

Accessory Code	Name
1903008011	Printer
1903008012	Printing set
1903008013	Load thermometer
1903008014	Cooling fan
1903008015	Temperature controller
1903008016	Temperature sensor
1903008017	Pressure Sensor
1903008018	Heater
1903008019	Sealing ring

Accessory Code	Name
1903008020	Pressure gauge
1903008021	Stainless steel pail
1903008022	Relay
1903008023	Main board
1903008024	Display board
1903008025	Jiggle switch
1903008026	Handle
1903008027	Pressure switch
1903008028	Solenoid Valve



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)