

PORTABLE STEAM STERILIZER LTAUT35-10



PORTABLE STEAM STERILIZER LTAUT35-10

Used in Medical, Research Institutions, Laboratory, Drinking Water Plant, Food industry, Hospital, Agriculture, Hospital, Scientific Research..

Also known as Portable Sterilizer, Steam Autoclave Sterilizer, Portable Autoclave, Pressure Steam Sterilizer, Portable Pressure Steam Sterilizer, Dental Autoclave.

LTAUT35-10 PORTABLE STEAM STERILIZER

SPECIFICATIONS

Model	LTAUT35-10
Capacity	18 L
Function	Electric heated
Sterilization Temperature	126°C
Steam Exhausting Temperature	≤1°C
Working Pressure	0.14-0.16Mpa
Maximum Working Pressure	0.165 Mpa
Timer Range	0-99min
Chamber Dimension	280x260mm
Overall Dimension (LxWxH)	440x420x630 mm
Chamber Material	Stainless Steel
Weight	12/14 kg
Power	2 KW
Power Supply	220 V, 50 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com