

BENCHTOP AUTOCLAVE LTAUT35-5



BENCHTOP AUTOCLAVE LTAUT35-5

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Medical Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

LTAUT35-5 BENCHTOP AUTOCLAVE



SPECIFICATIONS

Model	LTAUT35-5
Capacity	8 L
Sterilization Temperature	121°C, 134°C
Sterilization Time	0-60 min
Times of Sterilization after filled up	3-4 times
Working Pressure	0.21 - 0.23 MPa
Pressure Sensor	Yes
Precision of Pressure Sensor	0.1 KPA
Precision of Temperature Sensor	0.1°C
Maximum Vacuum in Bar	0.82
Water supply system	Manual
Quantity of clean water tank	1
Quantity of waste water tank	1
Water drain pipe	1
Quantity of water pump	1
Water Tank	2 x 1.5 L Internal Stainless Steel water tanks
Quantity of steam generator	1
Quantity of vacuum pump	1
Type of vacuum pump	Double head
Water Filter	1
Vapour Filter	Yes
Temperature Sensor	Yes
Automatic protection when overpressure	Yes
Automatic protection when over temperature	Yes
Bacteria Filter	1

Lock system of door	Mechanic + Electrical lockers
Screen size	2.5 x 2 cm, 3.8x 2 cm
Life span	5 years
Hot air exhaust system	Copper condenser, fan
121°C Solid	Yes
121°C Porous	Yes
121°C Hollow	Yes
134°C Solid	Yes
134°C Porous	Yes
134°C Hollow	Yes
B&D Test	Yes
Helix Test	No
Vacuum Test	Yes
Clean Program	No
Prion Program	No
User Define	No
Chamber Dimension	500x1000 mm
Diameter of chamber	170 mm
Depth of chamber	320 mm
Inner Dimension	535 X 365 X 320 mm
Overall Dimension (LxWxH)	1000x600x1200 mm
Minimum size of table for installation	350 X 280 mm
Weight	34/40 kg
Power	1700 W
Power Supply	220 V, 60 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com