

PORTABLE STEAM STERILIZER LTAUT35-6



PORTABLE STEAM STERILIZER LTAUT35-6

Used in Medical, Research Institutions, Laboratory, Drinking Water Plant, Food industry, Hospital, Agriculture, Hospital, Scientific Research..

Also known as Portable Sterilizer, Steam Autoclave Sterilizer, Portable Autoclave, Pressure Steam Sterilizer, Portable Pressure Steam Sterilizer, Dental Autoclave.

LTAUT35-6 PORTABLE STEAM STERILIZER



SPECIFICATIONS

Model	LTAUT35-6
Capacity	18 L
Function	Electric or LPG heated.
Sterilization Temperature	126°C
Working Pressure	0.14-0.16 Mpa
Maximum Working Pressure	0.165 Mpa
Chamber Dimension	280x260 mm
Overall Dimension (LxWxH)	410x410x430 mm
Chamber Material	Stainless Steel
Weight	14/16 kg
Power	2 KW
Power Supply	220 V, 50 Hz

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com