

# PORTABLE STEAM STERILIZER LTAUT35-9



# PORTABLE STEAM STERILIZER LTAUT35-9

Used in Medical, Research Institutions, Laboratory, Drinking Water Plant, Food industry, Hospital, Agriculture, Hospital, Scientific Research..

Also known as Portable Sterilizer, Steam Autoclave Sterilizer, Portable Autoclave, Pressure Steam Sterilizer, Portable Pressure Steam Sterilizer, Dental Autoclave.

## LTAUT35-9 PORTABLE STEAM STERILIZER



### SPECIFICATIONS

Model	LTAUT35-9
Capacity	24 L
Function	Electric heated
Sterilization Temperature	126°C
Steam Exhausting Temperature	≤1°C
Working Pressure	0.14-0.16Mpa
Maximum Working Pressure	0.165 Mpa
Chamber Dimension	280x390mm
Overall Dimension (LxWxH)	410x410x630 mm
Chamber Material	Stainless Steel
Weight	15/17 kg
Power	2 KW
Power Supply	220 V, 50 Hz



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)