

# DENSITY METER BALANCE LTBAL10-1



# DENSITY METER BALANCE LTBAL10-1

Used in Density Determination, Industrial, Food Retailing, Laboratory, Research.

Also known as Density Scales, Laboratory Density Balances, Laboratory Density Scales.

## LTBAL10-1 DENSITY METER BALANCE

Density frame  
Electronic Balance  
RS232 interface



### SPECIFICATIONS

Model	LTBAL10-1
Capacity (g)	310 g
Resolution (g)	0.01 g
Min Weighing (g)	0.4 g
Display	LCD(white back light with black font)
Sensor	Load Cell
Pan Size	160x70 mm
Operation Temperature	5-25°C
Repeatability	±0.02 g
Liner	±0.03 g
Cal. Weight	200 g
Test Mode	Solid, Liquid
Case Material	Plastic Case with Aluminium density kit
Unit Mode	Weighing: g / ct / oz Density: D
Overall Dimension	139x100x78 mm
Packing Size	440x310x300 mmx2
Gross Weight	4 kg

### OPTIONAL ACCESSORIES

Accessory Code	Name
3800808008	RS232C & RS485 Interface
3800808009	Printer



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)