

# DENSITY METER BALANCE LTBAL10-2



# DENSITY METER BALANCE LTBAL10-2

Used in Density Determination, Industrial, Food Retailing, Laboratory, Research.

Also known as Density Scales, Laboratory Density Balances, Laboratory Density Scales.

## LTBAL10-2 DENSITY METER BALANCE

Density frame  
Electronic Balance  
RS232 interface



### SPECIFICATIONS

Model	LTBAL10-2
Capacity (g)	210 g
Resolution (g)	0.001 g
Min Weighing (g)	0.004 g
Display	LCD(white back light with black font)
Sensor	Load Cell
Pan Size	160x70 mm
Operation Temperature	5-25°C
Repeatability	±0.002 g
Liner	±0.003 g
Cal. Weight	200 g
Test Mode	Solid, Liquid, Powder, Absorbent, etc
Case Material	Aluminium Case with Aluminium density kit
Unit Mode	Weighing: g / ct / oz Density: D / % / K
Overall Dimension	139x100x78 mm
Packing Size	440x310x300 mmx2
Gross Weight	6 kg

### OPTIONAL ACCESSORIES

Accessory Code	Name
3800808008	RS232C & RS485 Interface
3800808009	Printer



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)