

# ANALYTICAL BALANCE LTBAL8-7



# ANALYTICAL BALANCE LTBAL8-7

Used in Totalization, Formulation, SQC, Industry, Laboratory, Research.

Also known as Lab Balances, Laboratory Analytical Balances, Analytical Scales, Laboratory Analytical Scales, Weighing Scales.

## LTBAL8-7 ANALYTICAL BALANCE

Electromagnetic balanced sensors

The multiple functions modes with tare, counting, unit conversion and weighing

External Calibration



## SPECIFICATIONS

Model	LTBAL8-7
Capacity (g)	200 g
Resolution (g)	0.0001 g
Min Weighing (g)	0.0004 g
Stable Time	≤3 sec
Display	LCD(white back light with black font)
Pan Size	Φ 80 mm
Operation Temperature	5-25°C
Repeatability	±0.0002 g
Liner	±0.0003 g
Cal. Weight	200 g
Overall Dimension	240x190x265 mm
Packing Size	Inside: 415x285x430 mm Outside: 475x315x485 mm
Gross Weight	10 kg

## OPTIONAL ACCESSORIES

Accessory Code	Name
3800808008	RS232C & RS485 Interface
3800808009	Printer



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)