

TOPLOADING PRECISION BALANCE LTBAL9-2



TOPLOADING PRECISION BALANCE LTBAL9-2

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balances, Laboratory Precision Balances, Laboratory Loading Balances, Precision Scales.

LTBAL9-2 TOPLOADING PRECISION BALANCE

Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing

4 conversion units: oz,ct,lb,g



SPECIFICATIONS

| | |
|-----------------------|--|
| Model | LTBAL9-2 |
| Capacity (g) | 1100 g |
| Resolution (g) | 0.1 g |
| Min Weighing (g) | 0.4 g |
| Stable Time | ≤2 sec |
| Pan Size | Φ 130 mm |
| Operation Temperature | 5-25°C |
| Repeatability | 0.2 g |
| Liner | 0.3 g |
| Cal. Weight | 500 g |
| Packing Size | Inside: 310x230x130 mm Outside: 660x470x320 mm |
| Gross Weight | 2 / 22 kg |

OPTIONAL ACCESSORIES

| Accessory Code | Name |
|----------------|--------------------------|
| 3800808008 | RS232C & RS485 Interface |
| 3800808009 | Printer |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com