

# TOPLOADING PRECISION BALANCE LTBAL9-22



# TOPLOADING PRECISION BALANCE LTBAL9-22

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balances, Laboratory Precision Balances, Laboratory Loading Balances, Precision Scales.

## LTBAL9-22 TOPLOADING PRECISION BALANCE

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing



### SPECIFICATIONS

|                       |  |
|-----------------------|--|
| Model                 | LTBAL9-22                                      |
| Capacity (g)          | 310 g  |
| Resolution (g)        | 0.001 g  |
| Min Weighing (g)      | 0.004 g  |
| Stable Time           | ≤3 sec   |
| Display               | LCD(white back light with black font)          |
| Sensor                | Load Cell                                      |
| Pan Size              | Φ 90 mm  |
| Operation Temperature | 5-25°C   |
| Repeatability         | ±0.002 g                                       |
| Liner                 | ±0.003 g                                       |
| Cal. Weight           | 200 g  |
| Overall Dimension     | 177x173x205 mm                                 |
| Packing Size          | Inside: 350x265x345 mm Outside: 740x550x390 mm |
| Gross Weight          | 4 / 16 kg                                      |

### OPTIONAL ACCESSORIES

| Accessory Code | Name                     |
|----------------|--------------------------|
| 3800808008     | RS232C & RS485 Interface |
| 3800808009     | Printer                  |
| 3800927008     | USB                      |
| 3800927009     | Dual Display             |



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)