

# TOPLOADING PRECISION BALANCE LTBAL9-23



# TOPLOADING PRECISION BALANCE LTBAL9-23

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balances, Laboratory Precision Balances, Laboratory Loading Balances, Precision Scales.

## LTBAL9-23 TOPLOADING PRECISION BALANCE

Aluminium shell, Anti-jamming, Anti-radiation

LED digital display

The multiple functions modes with tare, counting, unit conversion and weighing



## SPECIFICATIONS

Model	LTBAL9-23
Capacity (g)	3100 g
Resolution (g)	0.01 g
Min Weighing (g)	0.04 g
Stable Time	3 sec
Display	LED / LCD
Sensor	Load Cell
Pan Size	180x180 mm
Operation Temperature	5-25°C
Repeatability	±0.03 g
Cal. Weight	2000 g
Gross Weight	5 kg

## OPTIONAL ACCESSORIES

Accessory Code	Name
3800808008	RS232C & RS485 Interface
3800808009	Printer
3800927008	USB



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)