

TOPLOADING PRECISION BALANCE LTBAL9-25



TOPLOADING PRECISION BALANCE LTBAL9-25

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balances, Laboratory Precision Balances, Laboratory Loading Balances, Precision Scales.

LTBAL9-25 TOPLOADING PRECISION BALANCE

Aluminium shell, Anti-jamming, Anti-radiation

LED digital display

The multiple functions modes with tare, counting, unit conversion and weighing



SPECIFICATIONS

Model	LTBAL9-25
Capacity (g)	5100 g
Resolution (g)	0.01 g
Min Weighing (g)	0.04 g
Stable Time	3 sec
Display	LED / LCD
Sensor	Load Cell
Pan Size	180x180 mm
Operation Temperature	5-25°C
Repeatability	±0.03 g
Cal. Weight	2000 g
Gross Weight	5 kg

OPTIONAL ACCESSORIES

Accessory Code	Name
3800808008	RS232C & RS485 Interface
3800808009	Printer
3800927008	USB



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com