

# TOPLOADING PRECISION BALANCE LTBAL9-3



# TOPLOADING PRECISION BALANCE LTBAL9-3

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balances, Laboratory Precision Balances, Laboratory Loading Balances, Precision Scales.

## LTBAL9-3 TOPLOADING PRECISION BALANCE

Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing

4 conversion units: oz,ct,lb,g



### SPECIFICATIONS

Model	LTBAL9-3
Capacity (g)	2100 g
Resolution (g)	0.1 g
Min Weighing (g)	0.4 g
Stable Time	≤2 sec
Pan Size	Φ 130 mm
Operation Temperature	5-25°C
Repeatability	0.2 g
Liner	0.3 g
Cal. Weight	1000 g
Packing Size	Inside: 310x230x130 mm Outside: 660x470x320 mm
Gross Weight	2 / 22 kg

### OPTIONAL ACCESSORIES

Accessory Code	Name
3800808008	RS232C & RS485 Interface
3800808009	Printer



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)