

# WATER BATH LTBAT10-1



# WATER BATH LTBAT10-1

Used in Clinical, General Chemistry, Preservation, Reactions of sample, DNA amplification, Initial Denaturation of electrophoresis, Laboratory.

Also known as Digital Heating Blocks, Laboratory Mini Dry Bath Incubators.

## LTBAT10-1 WATER BATH

Unique and novel appearance, small size, powerful functions. All surfaces in contact with water are made of high-quality stainless steel and are corrosion-resistant.

TFT LCD holographic displays operating speed, temperature, time. Graphic dynamic symbols and precise data at a glance. DC motor, low noise, maintenance-free, long service life.



## SPECIFICATIONS

Model	LTBAT10-1
Capacity	10 L
Temp. Control Range	R.T.+5°C ~100°C
Temp.Stability(100°C)	±0.5°C
Time Range	1m~99h59m OFF is forever
Oscillation mode	Reciprocating
Amplitude	26 mm
Speed range	20~200 rpm
Sink size((WxDxH)	250x304x200 mm
Voltage	AC220/110V, 50/60 Hz
Power	1100W
Dimensions(WxDxH)	W.288 x D.506 x H.350 mm
Net weight	10kgs



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)