

# WATER BATH LTBAT13-4



# WATER BATH LTBAT13-4

Used in Enzyme reactions, Immunoassay, Incubation, Laboratory, Research, Medical, Pharmaceuticals.

Also known as Laboratory water baths.

## LTBAT13-4 WATER BATH

Unique and novel appearance, small size, powerful functions.

All surface that contact with water adopt stainless steel material and PC material, with strong corrosion resistance.

With one hole and volume of 2 liter.

Using microcomputer program PID automatic temperature control, uniform temperature control, high precision.

Built-in water level sensor, to prevent the hidden dangers of dry burning, safe and reliable.



## SPECIFICATIONS

Model	LTBAT13-4
Time Range	1~99h59min (0 is forever)
Temperature Range	RT. +5°C~100°C
Temp. accuracy	±0.5°C
Temp. Uniformity@37°C	±0.5°C
Number of holes	3
Volume	2Lx3
Workspace	140x126x95mmx3
Voltage	AC 220V 50/60Hz
Power	900W
Dimensions	W.510xD.290xH.240mm
Net weight	6.7kgs



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)