

ELECTRIC HEATING THERMOSTATIC OIL BATH LTBAT21-1



ELECTRIC HEATING THERMOSTATIC OIL BATH

LTBAT21-1

Used in High temperature organic distillation, External circulation, Calibration tests, High temperature experiments, Laboratory.

Also known as Cooling Bath Circulators, Refrigerated Circulators, Laboratory Refrigerated Bath Circulators.

LTBAT21-1 ELECTRIC HEATING THERMOSTATIC OIL BATH

Advanced inner circulation and outer circulation pump

LED dual-window digital display

High effective thermal insulation materials



SPECIFICATIONS

Model	LTBAT21-1
Capacity	10 L
Temperature Range	RT +5 to 300°C
Temperature Accuracy	±0.1°C
Pump Flow Rate	8 L/min
Inner Dimension (WxDxH)	250x200x200 mm
Control	Microprocessor Controller
Display	LED
Power Supply	220 V

LABTARE

Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com