

DRY BATH INCUBATOR LTBAT21-14



DRY BATH INCUBATOR LTBAT21-14

Used in Calibration tests, High temperature organic distillation, External circulation, High temperature experiments.
Also known as Bath Circulators, Laboratory Bath Circulators.

LTBAT21-14 DRY BATH INCUBATOR

- Microprocessor temperature controller
- Easy maintenance
- LCD display with touch key panel for temperature setting
- Over temperature alarm system



SPECIFICATIONS

Model	LTBAT21-14
Temperature Range	RT+5-100°C
Display Accuracy	0.1°C
Temperature Control Accuracy	±0.3°C



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com