

DRY BATH INCUBATOR LTBAT21-15



DRY BATH INCUBATOR LTBAT21-15

Used in Calibration tests, High temperature organic distillation, External circulation, High temperature experiments.
Also known as Bath Circulators, Laboratory Bath Circulators.

LTBAT21-15 DRY BATH INCUBATOR

- Microprocessor temperature controller
- Easy maintenance
- LCD display with touch key panel for temperature setting
- Over temperature alarm system



SPECIFICATIONS

Model	LTBAT21-15
Temperature Range	Room temp.-120°C
Temperature Uniformity	±0.5°C
Outer Dimension (WxDxH)	95Wx150Dx50H mm
Weight	3 kg
Temperature Control Accuracy	±0.5°C
Overheating protection	140°C
Adapter block material	Aluminium
Display	LED
Power	160 W
Power Supply	100-120 V / 220-240 V ,50 Hz / 60 Hz



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com