

DRY BATH INCUBATOR LTBAT21-16



DRY BATH INCUBATOR LTBAT21-16

Used in Calibration tests, High temperature organic distillation, External circulation, High temperature experiments.
Also known as Bath Circulators, Laboratory Bath Circulators.

LTBAT21-16 DRY BATH INCUBATOR

Microprocessor temperature controller

Easy maintenance

LCD display with touch key panel for temperature setting

Over temperature alarm system



SPECIFICATIONS

Model	LTBAT21-16
Temperature Uniformity	$\pm 0.5^{\circ}\text{C}$
Outer Dimension (WxDxH)	185x185x90 mm
Weight	1.4 kg
Temperature Precision	$\pm 0.5^{\circ}\text{C}$
Operation Ambient Temperature	8-30°C
Power Supply	12 V DC, 5A



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com