

# AIR BATH LTBAT21-18



# AIR BATH LTBAT21-18

Used in High temperature organic distillation, External circulation, Calibration tests, High temperature experiments, Laboratory.

Also known as Cooling Bath Circulators, Refrigerated Circulators, Laboratory Refrigerated Bath Circulators.

## LTBAT21-18 AIR BATH



### SPECIFICATIONS

Model	LTBAT21-18
Temp. Control Range	R.T.+5°C~60°C
Temp. Setting Range	5°C~60°C
Temp. Stability@37°C	±0.2°C
Temp. Uniformity	±0.5°C
Temp. Display Accuracy	0.1°C
Time Range	1min~99h59min
Fan Speed Range	1200~3200rpm
Incubation Space Size	L.245 x W.215 x H.85mm
Voltage	AC 220V/AC 110V, 50/60Hz
Power	300W
Fuse	250 V, 2A/4A, φ5x20
Dimension	W.275 x D.395 x H.150mm
Net Weight	8.1 kgs



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)