

# ELECTRIC HEATING THERMOSTATIC OIL BATH LTBAT21-4



# ELECTRIC HEATING THERMOSTATIC OIL BATH

## LTBAT21-4

Used in High temperature organic distillation, External circulation, Calibration tests, High temperature experiments, Laboratory.

Also known as Cooling Bath Circulators, Refrigerated Circulators, Laboratory Refrigerated Bath Circulators.

## LTBAT21-4 ELECTRIC HEATING THERMOSTATIC OIL BATH

Advanced inner circulation and outer circulation pump

LED dual-window digital display

High effective thermal insulation materials



### SPECIFICATIONS

Model	LTBAT21-4
Capacity	50 L
Temperature Range	RT +5 to 300°C
Temperature Accuracy	±0.1°C
Pump Flow Rate	8 L/min
Inner Dimension (WxDxH)	400x350x350 mm
Control	Microprocessor Controller
Display	LED
Power Supply	220 V

# LABTARE

**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)