

DRY BATH INCUBATOR LTBAT21-8



DRY BATH INCUBATOR LTBAT21-8

Used in Pharmaceutical, Chemical, Food Safety, Environment, Quality inspection, Laboratory.
Also known as Laboratory Dual Temperature Dry Bath Incubators.

LTBAT21-8 DRY BATH INCUBATOR



SPECIFICATIONS

Model	LTBAT21-8
Capacity	10 L
Function	Power adjustable
Temperature Range	R.T.+5 °C~100 °C
Temperature Control Accuracy	±0.5°C
Cleaning Tank (LxWxH)	300x240x150 mm
Outer Dimension (WxDxH)	420x310x370 mm
Weight	3 kg
Maximum Input Power	400
Maximum Output Power	360
Heating Power	NO
Power Adjustable (%)	30-100
Drain	YES
Level Control System	YES



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com