

DRY BATH INCUBATOR LTBAT21-9



DRY BATH INCUBATOR LTBAT21-9

Used in Pharmaceutical, Chemical, Food Safety, Environment, Quality inspection, Laboratory.
Also known as Laboratory Dual Temperature Dry Bath Incubators.

LTBAT21-9 DRY BATH INCUBATOR



SPECIFICATIONS

Model	LTBAT21-9
Capacity	10 L
Function	Basic model
Temperature Range	RT+5-100 °C
Maximum Temperature Difference at 40°C	±0.5°C at 100°C
Display Accuracy	±0.1°C
Cleaning Tank (LxWxH)	300x240x150 mm
Outer Dimension (WxDxH)	420x310x370 mm
Weight	2.5 kg
Maximum Input Power	400
Maximum Output Power	360
Heating Power	NO
Power Adjustable (%)	NO
Drain	YES
Level Control System	YES



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com