

# THERMOSTATIC OIL BATH LTBAT22-1



# THERMOSTATIC OIL BATH LTBAT22-1

Used in Immunoassays, Melting points, Enzyme reactions, Boiling points, Incubation, Activation of cultures, Laboratory procedures..

Also known as Digital Dry Bath Incubators, Heating blocks, Laboratory Dry Bath Incubators.

## LTBAT22-1 THERMOSTATIC OIL BATH



### SPECIFICATIONS

Model	LTBAT22-1
Temperature Range	RT~100°C
Temperature Fluctuation	±0.1°C
Pump Flow Rate	8 L / min
Inner Dimension (WxDxH)	200x250x200 mm
Power Supply	220 V



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)