

THERMOSTATIC OIL BATH LTBAT22-2



THERMOSTATIC OIL BATH LTBAT22-2

Used in Immunoassays, Melting points, Enzyme reactions, Boiling points, Incubation, Activation of cultures, Laboratory procedures..

Also known as Digital Dry Bath Incubators, Heating blocks, Laboratory Dry Bath Incubators.

LTBAT22-2 THERMOSTATIC OIL BATH

- Microprocessor controlled incubation temperature and time
- Compact design with stable operation
- Rapid heat up, uniform heating, high stability, low energy consumption and noise
- Simultaneous display of set temperature and time
- Built in temperature calibration function
- Aluminium blocks ensures even heat distribution, eliminating the possibility of heater burnout
- Custom blocks available to meet experimental requirements
- Easy cleaning, replacement and disinfection of metal blocks
- Automatic fault detection and buzzer alarm function
- Audio alarm indicates program completion
- Over temperature protection device



SPECIFICATIONS

| | |
|-------------------------|----------------|
| Model | LTBAT22-2 |
| Temperature Range | RT~100 °C |
| Temperature Fluctuation | ±0.1°C |
| Pump Flow Rate | 8 L / min |
| Inner Dimension (WxDxH) | 280x400x180 mm |
| Power Supply | 220 V |



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com