

THERMOSTATIC OIL BATH LTBAT22-3



THERMOSTATIC OIL BATH LTBAT22-3

Used in Immunoassays, Melting points, Enzyme reactions, Boiling points, Incubation, Activation of cultures, Laboratory procedures..

Also known as Digital Dry Bath Incubators, Heating blocks, Laboratory Dry Bath Incubators.

LTBAT22-3 THERMOSTATIC OIL BATH

Microprocessor controlled incubation temperature and time

Compact design with stable operation

Rapid heat up, uniform heating, high stability, low energy consumption and noise

Simultaneous display of set temperature and time

Built in temperature calibration function

Aluminium blocks ensures even heat distribution, eliminating the possibility of heater burnout

Custom blocks available to meet experimental requirements

Easy cleaning, replacement and disinfection of metal blocks

Automatic fault detection and buzzer alarm function

Audio alarm indicates program completion

Over temperature protection device



SPECIFICATIONS

Model	LTBAT22-3
Temperature Range	RT~100 °C
Temperature Fluctuation	±0.1°C
Pump Flow Rate	8 L / min
Inner Dimension (WxDxH)	300x400x250 mm
Power Supply	220 V



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com