

DRY BATH INCUBATOR LTBAT23-19



DRY BATH INCUBATOR LTBAT23-19

Used in Clinical, General Chemistry, Preservation, Reactions of sample, DNA amplification, Initial Denaturation of electrophoresis, Laboratory.

Also known as Digital Heating Blocks, Laboratory Mini Dry Bath Incubators.

LTBAT23-19 DRY BATH INCUBATOR



SPECIFICATIONS

Model	LTBAT23-19
Temperature Range	-10°C-100°C
Temperature Accuracy for lid	±1°C
Maximum Temperature Difference at 40°C	±0.5°C at 100°C
Display Accuracy	±0.1°C
Outer Dimension (WxDxH)	300Lx220Wx180H mm
Weight	5 kg



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com