

DRY BATH INCUBATOR LTBAT23-20



DRY BATH INCUBATOR LTBAT23-20

Used in Clinical, General Chemistry, Preservation, Reactions of sample, DNA amplification, Initial Denaturation of electrophoresis, Laboratory.

Also known as Digital Heating Blocks, Laboratory Mini Dry Bath Incubators.

LTBAT23-20 DRY BATH INCUBATOR



SPECIFICATIONS

Model	LTBAT23-20
Temperature Range	0°C-100°C if Ambient Temperature ≤ 20°C 4°C-100°C if Ambient Temperature ≤ 25°C 10°C-100°C if Ambient Temperature ≤ 30°C
Temperature Accuracy for lid	±1°C
Display Accuracy	0.1°C
Temperature Control Accuracy	≤ ±0.3°C



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com