

# MINI DRY BATH INCUBATOR LTBAT24-10



# MINI DRY BATH INCUBATOR LTBAT24-10

Used in High temperature organic distillation, External circulation, Calibration tests, High temperature experiments, Laboratory.

Also known as Cooling Bath Circulators, Refrigerated Circulators, Laboratory Refrigerated Bath Circulators.

## LTBAT24-10 MINI DRY BATH INCUBATOR

Fuzzy PID temperature control system: control temperature stability accurately and fast.

Multi-protection: Over-current protection, high and low temperature warning functions, dry-running protection.

Convenience function: Parameters memory function, temperature correction.

ECO-refrigerant R134a and high quality compressor guarantee the cooling system ecofriendly and powerful

Low Noise design: use high quality compressor and mute pump.

Robust water tank: ultra-slushing and rust-proof

Bright LCD temperature display

Compact design, removable venting grid for convenient cleaning to maintain cooling performance and drain off water



## SPECIFICATIONS

|                                     |                   |
|-------------------------------------|-------------------|
| Model                               | LTBAT24-10        |
| Temperature Range                   | RT+5-100°C        |
| Display Accuracy                    | ±0.1°C            |
| Outer Dimension (WxDxH)             | 152Wx120Dx112H mm |
| Weight                              | 0.9 kg            |
| Lid's Temperature range             | RT+5-100°C        |
| Temperature Control Accuracy        | ±0.5°C            |
| Temperature Control Accuracy @ 40°C | ±0.3°C            |
| Power Supply                        | 24 V DC           |



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)