

MINI DRY BATH INCUBATOR LTBAT25-13



MINI DRY BATH INCUBATOR LTBAT25-13

Used in Clinical, General Chemistry, Preservation, Reactions of sample, DNA amplification, Initial Denaturation of electrophoresis, Laboratory.

Also known as Digital Heating Blocks, Laboratory Mini Dry Bath Incubators.

LTBAT25-13 MINI DRY BATH INCUBATOR



SPECIFICATIONS

Model	LTBAT25-13
Temperature Range	RT+5-99.9°C (BX7103 & BX7104) 0-99.9°C(BX7102)
Temperature Stability	0.2°C
Display Accuracy	0.1°C
Outer Dimension (WxDxH)	360x180x150 mm
Weight	3 kg
Temperature Control Accuracy	≤0.5°C
Power	120 W
Power Supply	12 V DC
Note	* only for BX7102



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com