

# MINI DRY BATH LTBAT25-3



# MINI DRY BATH LTBAT25-3

Used in Immunoassays, Melting points, Enzyme reactions, Boiling points, Incubation, Activation of cultures, Laboratory procedures..

Also known as Digital Dry Bath Incubators, Heating blocks, Laboratory Dry Bath Incubators.

## LTBAT25-3 MINI DRY BATH



### SPECIFICATIONS

Model	LTBAT25-3
Capacity	4 L
Function	Heating function
Temperature Range	+5°C~100°C
Temperature Uniformity	≤±0.3°C
Cleaning Tank (LxWxH)	300x150x100 mm
Outer Dimension (WxDxH)	420x230x270 mm
Weight	0.5 kg
Maximum Input Power	450
Maximum Output Power	180
Heating Power	250
Power Adjustable (%)	30-100
Drain	NO
Level Control System	NO
Power Supply	AC 100-240 V, 50 / 60 Hz



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)