

MINI DRY BATH INCUBATOR LTBAT25-4



MINI DRY BATH INCUBATOR LTBAT25-4

Used in Immunoassays, Melting points, Enzyme reactions, Boiling points, Incubation, Activation of cultures, Laboratory procedures..

Also known as Digital Dry Bath Incubators, Heating blocks, Laboratory Dry Bath Incubators.

LTBAT25-4 MINI DRY BATH INCUBATOR



SPECIFICATIONS

Model	LTBAT25-4
Temperature Range	RT+5-80 °C
Display Accuracy	±0.1 °C
Weight	0.5 kg
Temperature Control Accuracy	±0.5°C
Temperature Control Accuracy @ 40°C	±0.3°C
Power Supply	12 V DC



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com