

MINI DRY BATH INCUBATOR LTBAT26-6



MINI DRY BATH INCUBATOR LTBAT26-6

Used in Immunoassays, Melting points, Enzyme reactions, Boiling points, Incubation, Activation of cultures, Laboratory procedures..

Also known as Digital Dry Bath Incubators, Heating blocks, Laboratory Dry Bath Incubators.

LTBAT26-6 MINI DRY BATH INCUBATOR



SPECIFICATIONS

Model	LTBAT26-6
Temperature Uniformity	$\pm 0.5^{\circ}\text{C}$
Weight	0.75 kg
Temperature Precision	$\pm 0.5^{\circ}\text{C}$
Operation Ambient Temperature	8-30°C
Power Supply	12 V DC, 5A



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com