

DRY BATH CHILLER LTBAT26-7



DRY BATH CHILLER LTBAT26-7

Used in Pharmaceutical, Chemical, Food Safety, Environment, Quality inspection, Laboratory.
Also known as Laboratory Dual Temperature Dry Bath Incubators.

LTBAT26-7 DRY BATH CHILLER



SPECIFICATIONS

Model	LTBAT26-7
Capacity	4 L
Function	Basic model
Temperature Range	Maintains at 0°C
Temperature Uniformity	±0.5°C
Cleaning Tank (LxWxH)	300x150x100 mm
Outer Dimension (WxDxH)	420x230x270 mm
Weight	1.4 kg
Temperature Precision	±0.5°C
Operation Ambient Temperature	8-30°C
Maximum Input Power	220
Maximum Output Power	180
Heating Power	NO
Power Adjustable (%)	NO
Drain	NO
Level Control System	NO
Power Supply	12 V DC, 7A



Labtare Analytical Instruments

Email: info@labtare.com | Website: labtare.com