

# LOW TEMPERATURE CIRCULATOR LTBAT30-2



# LOW TEMPERATURE CIRCULATOR LTBAT30-2

## LTBAT30-2 LOW TEMPERATURE CIRCULATOR

Microprocessor temperature controller with LCD screen ensures precise and reliable control, easy to operate

No angle in bath chamber, easy to clean

With interface to external water bath

Easy to set adjustable timer. (1 minute to 5,999 minutes)

Maintenance-free operation with easy to clean surface

Both heater and bath chamber are made of corrosion-resistant stainless steel

Drain valve makes emptying of bath water fast and easy for cleaning and moving

Programming setting function with 7 periods and 9 steps for each period, which means there are 63 programmable steps in total



### SPECIFICATIONS

Model	LTBAT30-2
Capacity	10.9 L
Temperature range	RT+5 - 150°C
Precision	±0.1
Liquid tank opening depth	240x170/150 mm
Power requirements	AC 220 V 50 HZ
Pump	8 L/min
Power Consumption	1050 W



**Labtare Analytical Instruments**

Email: [info@labtare.com](mailto:info@labtare.com) | Website: [labtare.com](http://labtare.com)